

KnowItAll Software Training

GC-MS Analysis Using KnowItAll GC Expert And SearchIt

GC-MS Analysis By Automatic Deconvolution

How to Perform GC-MS Analysis by Automatic Deconvolution

Purpose

These exercises demonstrate how to use KnowItAll GC Expert to automatically deconvolute/analyze GC-MS.

Objectives

These exercises will teach you:

- How to use KnowItAll GC Expert to auto-deconvolute GC-MS data into chemical component MS spectra which are automatically searched against millions of references.
 - How to generate reports.
-

Background

GC-MS data are information rich. Analysis can be time consuming, especially when examining complex analytes. We present a computer system that combines fast, flexible automated deconvolution and automatic database search to identify knowns and unknowns. Novel compounds can be identified and structural characteristics can be deduced from applying the MS Adaptive search that uses fragmentation and structural data to propose likely structural details of the unknown.

Training Files Used in This Lesson

- C:\Users\Public\Documents\Wiley\KnowItAll\Samples\GC Expert folder files

KnowItAll Applications Used

- KnowItAll GC Expert

GC-MS Deconvolution Algorithms

Our system follows individual m/z values across multiple spectra and extracts a pure spectrum from the data for each individual component while trying to separate components with overlapping m/z value peaks. If accurate m/z value data are available and the user selects this data to use it instead of unit m/z values, the chosen instrument accuracy (automatic, ppm, or fixed value) is used to determine the correct accurate m/z values throughout the entire GC-MS analysis. The m/z values in the raw data are converted into correct accurate values based on the closest value found by taking the instrument resolution into account. The corrected m/z values form the basis of the following deconvolution.

During the deconvolution step, individual m/z values are traced across multiple raw spectra, and a component spectrum is extracted while trying to separate components with overlapping m/z value peaks. The details of the algorithm are summarized to a large extent by the papers listed below¹⁻⁴.

Additional steps are added to automatically detect components with low intensity Reconstructed Total Ion Current (RTIC) chromatographic peaks as long as they can well be separated from neighboring components.

Although the algorithm is very complex, the details of the algorithm are summarized to a large extent by the following papers:

1. S. E. Stein. An Integrated Method for Spectrum Extraction and Compound Identification from Gas Chromatography/Mass Spectrometry Data. *J Am Soc Mass Spectrom* 1999, **10**, 770 –781.
2. R. G. Dromey, M. J. Stefik, T. C. Rindfleisch, A. M. Duffield. Extraction of Mass Spectra Free of Background and Neighboring Component Contributions from Gas Chromatography /Mass Spectrometry Data. *Analytical Chemistry*, 1976, **48(9)**, 1368-1375.
3. J. E. Biller, K. Biemann. Reconstructed Mass Spectra, A Novel Approach For The Utilization Of Gas Chromatography-Mass Spectrometer Data. *Analytical Letters* 1974, **7**, 515-28.
4. B. N. Colby. Spectral Deconvolution for Overlapping GC/MS Components. *J Am Soc Mass Spectrom* 1992, **3**, 558-562.

MS Spectral Comparison Algorithms

Research article

Journal of
MASS
SPECTROMETRY

Received: 5 October 2014

Revised: 16 February 2015

Accepted: 5 March 2015

Published online in Wiley Online Library

(wileyonlinelibrary.com) DOI 10.1002/jms.3591

Evaluation of mass spectral library search algorithms implemented in commercial software

Andrey Samokhin,^{a*} Ksenia Sotnezova,^a Vitaly Lashin^b and Igor Revelsky^a

MS SEARCH

Composite algorithm

$$SI = \frac{N_U \cdot \left[\frac{\left(\sum W_L \cdot W_U \right)^2}{\sum W_L^2 \sum W_U^2} \right] + \left[\sum \left(\frac{R_U}{R_L} \right)^n \right]}{N_U + N_{U\&L}}$$

Spectrum search type – identity (normal)

Presearch – default

Included Libs – MainLib

Apply limits – unchecked

Use constraints – unchecked

Dot-product algorithm

$$SI = \frac{\left(\sum W_L \cdot W_U \right)^2}{\sum W_L^2 \sum W_U^2}$$

Spectrum search type – similarity (simple)


Presearch – default

Included Libs – MainLib

Apply limits – unchecked

Use constraints – unchecked

Samokhin, K. Sotnezova, V. Lashin, I. Revelsky. Evaluation of mass spectral library search algorithms implemented in commercial software. *J. Mass Spectrom.* 2015, **50**, 820-825.



KnowItAll has 4 different algorithms:

Where

- Dot-Product (Cosine) – second equation in above image
 - Wiley Dot-Product (Cosine) (old KnowItAll algorithm) - the old Finnigan algorithm that verified at least 12 of the largest 16 peaks AND the base peak match before continuing with the dot product calculation.
 - Composite P1 – first equation in above image
 - Composite P3 – first equation in above image
- P1 and P3 are different by the power applied to the weighted intensity of the peak.

Example 1: GC-MS of Unit m/z Values

GUI explanation

The image below shows the deconvoluted GC-MS data of unit m/z values, and the search results for each component within the massive Wiley GC-MS database.

A – GC panel, showing the deconvoluted component peaks. One can use the checkboxes on the right side of the GC panel to turn on/off chromatogram(s) display.

When the box corresponding to the selected component peak of the chromatogram is selected, the box has a darker color;


B – Component Profile panel, highlights selected ions or components of the chromatogram(s);

C – Component Models panel, identifies the reference ion(s) used to model a component;

D – Database Match panel, shows the extracted spectrum (top) vs reference spectrum (bottom);

E – Components panel, where the Score is the combined spectrum search and reverse search Hit Quality Index (HQI) and, each component's GC area under curve (AUC) value;

G – Panel where parameters can be adjusted in algorithm.



Features of GC Expert

The Component Table (E) contains useful information:

- It shows each Components' estimated area.
- One can add **Notes** to a match. This can be done by clicking on the space below **Notes**, in the row of the matched component of interest.
- Right mouse-click brings up actions that the user can take. Among many, **Edit Component Columns** allows one to rearrange the display.

Components									
RT [min]	#	Match	Score	HQI	R.H...	Notes	Area %	Area	CAS Registry Number
15.9260	1	Bicyclo[4.3.1]decan-8-one, 1-bromo-	92.34	91.88	96.52	Interesting	1.40	1675158	56-23-5
16.5463	1	Ethane, 2-bromo-1,1-dichloro-	80.98	79.86	91.07		1.41	1685546	75-27-4
17.4897	1		67.37	65.48	84.40		0.41	491611	563-80-4
17.9705	1	Propanoic acid, 2,2-dichloro-	31.89	31.28	37.37		2.14	2558302	78-87-5
18.2397	1	1,3-Propanedithiol	35.82	34.32	49.30		1.69	2020873	10061-01-5
18.7663	1	3-(Chloromethyl)thiophene	30.20	28.74	43.33		1.95	2323999	79-01-6
19.2624	1	(3E)-1,3,7-Octatrien-5-yne	37.90	35.83	56.47		2.92	3490076	329-71-5
19.5835	1	1,4-Bis(ethenyl)benzene	51.34	50.52	58.70		2.91	3475740	124-48-1
19.6127	1	4-Chloroaniline						3621261	124-48-1
19.6533	1	(Trimethylsilyl)acetylene						2365877	79-00-5
20.7145	1	Propanoic acid, 2,2-dichloro-						2490372	110-75-8
21.0187	1	Myristic acid, methyl ester						252045	623-42-7
21.9853	1	3-Pentenoic acid, 4-methyl-						2344254	540-36-3
22.5330		No match found						1333680	
22.8777	1	9H-Fluoren-9-one						1703804	108-10-1
23.7254	1	2-Octanol, acetate						127023	2051-50-5
24.5232	1	Acetone-oxime						1652274	
24.8421		No match found						3787639	
24.8642		No match found						3078006	
24.9644	1	Dichloromethyl ethyl sulfone						1572112	
		<div><div>Remove Selected Component</div><div>Rank Selected Match Highest</div><div>View Component Hit List in Minelt</div><div>Edit Component Table Columns...</div><div>Copy Component Information</div><div>Copy All Component Information</div><div>Floating</div><div><input checked="" type="checkbox"/> Docking</div><div>Tabbed Document</div><div>Auto Hide</div></div>							
Component Models		Hide							

A hidden feature is the **Lookup** of a displayed property by using the **Control + F** action. For example, one can use **Lookup** to find a match to CAS Numbers such as CAS Number 329-7105.

Components									
RT [min]	#	Match	Score	HQI	R.H...	Notes	Area %	Area	CAS Registry Number
15.9260	1	Bicyclo[4.3.1]decan-8-one, 1-bromo-	92.34	91.88	96.52	Interesting	1.40	1675158	56-23-5
16.5463	1	Ethane, 2-bromo-1,1-dichloro-	80.98	79.86	91.07		1.41	1685546	75-27-4
17.4897	1		67.37	65.48	84.40		0.41	491611	563-80-4
17.9705	1	Propanoic acid, 2,2-dichloro-	31.89	31.28	37.37		2.14	2558302	78-87-5
18.2397	1	1,3-Propanedithiol	35.82	34.32	49.30		1.69	2020873	10061-01-5
18.7663	1	3-(Chloromethyl)thiophene	30.20	28.74	43.33		1.95	2323999	79-01-6
19.2624	1	(3E)-1,3,7-Octatrien-5-yne	37.90	35.83	56.47		2.92	3490076	329-71-5
19.5835	1	1,4-Bis(ethenyl)benzene	51.34	50.52	58.70		2.91	3475740	124-48-1
19.6127	1	4-Chloroaniline	35.77	34.29	49.11		3.03	3621261	124-48-1
19.6533	1	(Trimethylsilyl)acetylene	40.12	37.81	60.96		1.98	2365877	79-00-5
20.7145	1	Propanoic acid, 2,2-dichloro-	43.45	43.17	46.02		2.09	2490372	110-75-8
21.0187	1						0.21	252045	623-42-7
21.9853	1						1.96	2344254	540-36-3
22.5330							1.12	1333680	
22.8777	1						1.43	1703804	108-10-1
23.7254	1						0.11	127023	2051-50-5
24.5232	1						1.38	1652274	
24.8421							3.17	3787639	
24.8642							2.58	3078006	
24.9644	1	Dichloromethyl ethyl sulfone	70.37	67.89	94.71		1.32	1572112	

Find

Find what: 329-71-5



Find Next



Direction



☐ Match case
 ☐ Up
 ☒ Down



Cancel











Now, we will go through a complete process.



	Action	Result
1	<p>Go to the GC Expert application under Spectral Processing toolbar</p>  <p>Click Open Raw GC-MS Data File button</p> <p>Navigate to folder C:\Users\Public\Documents\Wiley\KnowItAll\Samp es\GC Expert\V50 - James Little.cdf file</p>	<p>GC Expert automatically performs:</p> <ol style="list-style-type: none"> 1. GC deconvolution to component extracted MS 2. Search extracted MS again reference databases 3. Report top hits 


	Action	Result
2	<p>Examine the hit lists for each deconvoluted components in the Components table</p> <p>Go to RT (MIN) 19.6127, notice the Score for this hit is very low, 56.74.</p>	 <p>From the comparison of Extracted Spectrum (top) and reference spectrum of methane, dibromochloro- in the bottom pane, we can see that the Extracted Spectrum contains methane, dibromochloro- as well as other components.</p>
3	<p>Transfer the Extracted Spectrum to SearchIt</p> <p>New in 2026 release: Reference databases from GC Expert are sent to SearchIt as the default databases to use when a selected component MS spectrum is transferred from GC Expert to SearchIt.</p>	

	Action	Result
4	<p>Only check the Extracted Spectrum at prompt.</p> <p>Click OK.</p>	
5	<p>In SearchIt, set to search for 2 (Mixture).</p>	

	Action	Result
6	<p>Note: New searching features in 2026 release.</p>	<p>Retention Index: one can optionally put RI (Retention Index) constraints:</p> <ul style="list-style-type: none"> • % Tolerance – the % difference of two component's Retention Index values • Low and High values: Retention Index value range  <p>The screenshot shows the SearchIt interface. Under 'Search Categories', 'Spectrum MS (GC)' is checked. Under 'RI Constraints', 'RI Constraints' is checked. The 'Property' dropdown is set to 'Retention Index'. The 'Tolerance' is set to 5%. The 'Low' and 'High' fields are empty.</p> <p>Partial Search: the spectrum search will not penalize peaks that are missing from the query compared to the database spectrum in the final HQI ("P.HQI").</p> <ul style="list-style-type: none"> • Select the checkbox below in SearchIt to execute. • Observe partial spectrum search results in MinIt as new P.HQI column. • Partial Search and Reverse Search cannot be used together.  <p>The screenshot shows the bottom section of the SearchIt interface. 'Number of components' is set to 1 (Single). 'Search Method' is set to Dot-Product (Cosine). 'Optimized Corrections' is checked. 'Accurate Mass Search' and 'Adaptive Search' are unchecked. 'Molecular m/z' is empty. 'Reverse Search' is unchecked. 'Partial Search' is checked and highlighted with a red box.</p>

	Action	Result																																											
7	Search.	<div>Mixture analysis returns a very good Composite Spectrum match top of the hit list:</div> <div><p>Below the plot is a table of search results:</p><table><tr><th>Table</th><th>Plot</th><th>Related Compounds View</th></tr><tr><th>HQI</th><th>Ratio</th><th>Exclude</th><th>Cor DB</th><th>ID</th><th>Name</th><th>Chemical Structure</th><th>Spectrum</th></tr><tr><td>1</td><td>94.09</td><td>N.A.</td><td></td><td></td><td>Composite Spectrum</td><td></td><td></td></tr><tr><td></td><td>0.59</td><td></td><td></td><td>WMSR3X</td><td>4571</td><td>1-Propene, 1,3-dichloro-, (Z)-</td><td></td></tr><tr><td></td><td>0.41</td><td></td><td></td><td>MSRX</td><td>27959</td><td>Methane, dibromochloro-</td><td></td></tr><tr><td></td><td>N.A.</td><td></td><td></td><td></td><td>Residual Spectrum</td><td></td><td></td></tr></table></div> <div>The Composite Spectrum in the first row is a mix of two components: 1-propene, 1,3-dichloro- (new component) and methane, dibromochloro-. The Residual Spectrum in the last row, the difference between Extracted Spectrum and Composite Spectrum is negligible, indicating no more components remain undetected.</div>	Table	Plot	Related Compounds View	HQI	Ratio	Exclude	Cor DB	ID	Name	Chemical Structure	Spectrum	1	94.09	N.A.			Composite Spectrum				0.59			WMSR3X	4571	1-Propene, 1,3-dichloro-, (Z)-			0.41			MSRX	27959	Methane, dibromochloro-			N.A.				Residual Spectrum		
Table	Plot	Related Compounds View																																											
HQI	Ratio	Exclude	Cor DB	ID	Name	Chemical Structure	Spectrum																																						
1	94.09	N.A.			Composite Spectrum																																								
	0.59			WMSR3X	4571	1-Propene, 1,3-dichloro-, (Z)-																																							
	0.41			MSRX	27959	Methane, dibromochloro-																																							
	N.A.				Residual Spectrum																																								


	Action	Result
8	<p>Go back to GC Expert,</p> <p>One can go to File > Settings menu to adjust component hit list parameters.</p>	<p>Hits evaluation parameters:</p>  <p>Reference database choices</p> 

	Action	Result
9	<p>One can click the TIC pane to add peaks manually, then select Analysis > Add Selected Peak To Analysis to perform a new analysis</p> <p>One can select View Component Hit List in Minelt to see Component table contents in Minelt</p>	 <p>The screenshot displays the KnowItAll software interface. The top menu bar includes File, Edit, View, Analysis, License, and Help. The Analysis menu is open, showing options: Masses, Add Selected Peak to Analysis, Set Baseline For Manually Selected Component, Use Manually Selected Slice As Component, Remove Selected Component, Rank Selected Match Highest, View Component Hit List in Minelt, Re-run Calculations, and Transfer to NIST MS Search... (Ctrl+Alt+N). On the left, a sidebar shows 'Spectral Processing' with sub-items 'Basics', 'Data', and 'Spectral Process'. Below this is a 'ProcessIt' button with a red 'X' icon.</p>

Example 2: High-Resolution GC-MS

GUI explanation




The instrument resolving power is required for the algorithm to automatically calculate the accurate m/z value of the data. Therefore, we use what our research considers a reasonable value by default that has a constant value and a variable value depending on mass (ppm). Empirically, this works in most cases. Increasing the m/z value accuracy too much incurs the danger of splitting an individual m/z value into individual mass spectral peaks that should be considered as just one. Decreasing the m/z value accuracy too much may cause individual mass spectral peaks to merge, resulting in incorrectly reported accurate m/z values.






If a user knows the instrument's resolving power, that value should be entered in the **Instrument Resolution** panel, highlighted by the purple box, in the above figure by clicking on 'Advanced' and manually entering the value.


It is possible to save the instrument's resolving power as part of a profile, which can then be selected depending on the type of instrument. A user can create several profiles with different resolution settings for different types of data (and instruments).

Hi-Resolution GC-MS Analysis using GC Expert – Exercise 1

	Action	Result
1	<p>Note: If you have changed your File > Settings in the previous exercise, you should go back to File > Settings to reset for this exercise.</p> <p>Start a new analysis by clicking on the Open Raw GC-MS Data File button.</p>  <p>Navigate to the folder: C:\Users\Public\Documents\Wiley\KnowItAll\Samples\GC Expert</p> <p>Select the file: Centroid Orbitrap HiRez - James Little.cdf</p> <p><i>Note: A pop-up window will appear. Click 'OK' to ignore the warning.</i></p> <p>Using the Components panel, find the RT (MIN) 6.9482 and click on that component.</p>	  <p><i>The m/z value of 170.0717 is slightly off for this compound.</i></p>

	Action	Result
2	<p>Check Advanced in the Instrument Resolution panel.</p> <p>Set the instrument resolution at 3 ppm,</p>  <p>Click OK.</p>	<p>GC Expert recognizes that this is a high-resolution data file and uniquely allows the user to fill in resolution information or adjust the resolution:</p> 
3		<p>The m/z value is corrected to 170.0726. This is the correct value. This instrument resolution set-up is persistent.</p> 

	Action	Result																																																																																																																																																																								
4	<p>To expand the component hits of a selected match, click on the (⊕) to the left of the component name.</p> <p>For example, find the component with RT 7.4282. Click the (⊕) icon.</p>	<div><div>Components</div><table><tr><th>RT [min]</th><th>#</th><th>Match</th><th>Score</th><th>HQI</th><th>R.H...</th><th>Com...</th><th>N</th></tr><tr><td>3.8140</td><td>1</td><td>Trimethyl-phosphate</td><td>89.17</td><td>89.12</td><td>89.61</td><td></td><td></td></tr><tr><td>3.8952</td><td>1</td><td>alpha-Pinene</td><td>91.44</td><td>91.41</td><td>91.73</td><td></td><td></td></tr><tr><td>6.8448</td><td>1</td><td>Biphenyl</td><td>97.09</td><td>96.96</td><td>98.27</td><td></td><td></td></tr><tr><td>6.9482</td><td>1</td><td>C-Phenyl-C-formyl-N-(2-furyl)nitron</td><td>90.20</td><td>89.52</td><td>96.24</td><td></td><td></td></tr><tr><td>6.9492</td><td>1</td><td>Diphenylether</td><td>90.85</td><td>90.77</td><td>91.50</td><td></td><td></td></tr><tr><td>7.0295</td><td>1</td><td>Diphenylether</td><td>76.30</td><td>75.40</td><td>84.41</td><td></td><td></td></tr><tr><td>7.0370</td><td>1</td><td>1,1'-Biphenyl, 3-ethoxy-</td><td>68.56</td><td>68.19</td><td>71.88</td><td></td><td></td></tr><tr><td>7.1110</td><td>1</td><td>Biphenyl</td><td>55.72</td><td>51.79</td><td>91.15</td><td></td><td></td></tr><tr><td>7.2828</td><td>1</td><td>Phenol, 3-(1,1-dimethylethyl)-4-metho...</td><td>93.02</td><td>92.97</td><td>93.49</td><td></td><td></td></tr><tr><td>7.3261</td><td>1</td><td>3-tert-Butyl-4-hydroxyanisol</td><td>95.83</td><td>95.80</td><td>96.11</td><td></td><td></td></tr><tr><td>7.4200</td><td>1</td><td>Dimethylisophthalate</td><td>90.00</td><td>89.93</td><td>90.65</td><td></td><td></td></tr><tr><td>7.4282</td><td>1</td><td>Butylated hydroxytoluene</td><td>93.35</td><td>93.25</td><td>94.26</td><td></td><td></td></tr><tr><td></td><td>2</td><td>Phenol, 2,6-bis(1,1-dimethylethyl)-4-...</td><td>90.14</td><td>90.05</td><td>90.96</td><td></td><td></td></tr><tr><td></td><td>3</td><td>Phenol, 2,6-bis(1,1-dimethylethyl)-4-...</td><td>90.09</td><td>89.98</td><td>91.07</td><td></td><td></td></tr><tr><td></td><td>4</td><td>Phenol, 2,4,6-tris(1-methylethyl)-</td><td>88.63</td><td>88.51</td><td>89.75</td><td></td><td></td></tr><tr><td></td><td>5</td><td>Butylated hydroxy toluene</td><td>87.63</td><td>87.09</td><td>92.45</td><td></td><td></td></tr><tr><td></td><td>6</td><td>Phenol, 2,4-bis(1,1-dimethylethyl)-6-...</td><td>86.72</td><td>86.64</td><td>87.44</td><td></td><td></td></tr><tr><td></td><td>7</td><td>(5S,6S)-6-(cyclohepta-2,4,6-trien-1-yl)-...</td><td>86.26</td><td>85.78</td><td>90.55</td><td></td><td></td></tr><tr><td></td><td>8</td><td>4-Hydroxy-3,5-diisopropylbenzaldehy...</td><td>81.40</td><td>81.34</td><td>81.93</td><td></td><td></td></tr><tr><td></td><td>9</td><td>Ethanone, 1-(5,6,7,8-tetrahydro-2,8,8-...</td><td>78.30</td><td>78.25</td><td>78.77</td><td></td><td></td></tr></table></div>	RT [min]	#	Match	Score	HQI	R.H...	Com...	N	3.8140	1	Trimethyl-phosphate	89.17	89.12	89.61			3.8952	1	alpha-Pinene	91.44	91.41	91.73			6.8448	1	Biphenyl	97.09	96.96	98.27			6.9482	1	C-Phenyl-C-formyl-N-(2-furyl)nitron	90.20	89.52	96.24			6.9492	1	Diphenylether	90.85	90.77	91.50			7.0295	1	Diphenylether	76.30	75.40	84.41			7.0370	1	1,1'-Biphenyl, 3-ethoxy-	68.56	68.19	71.88			7.1110	1	Biphenyl	55.72	51.79	91.15			7.2828	1	Phenol, 3-(1,1-dimethylethyl)-4-metho...	93.02	92.97	93.49			7.3261	1	3-tert-Butyl-4-hydroxyanisol	95.83	95.80	96.11			7.4200	1	Dimethylisophthalate	90.00	89.93	90.65			7.4282	1	Butylated hydroxytoluene	93.35	93.25	94.26				2	Phenol, 2,6-bis(1,1-dimethylethyl)-4-...	90.14	90.05	90.96				3	Phenol, 2,6-bis(1,1-dimethylethyl)-4-...	90.09	89.98	91.07				4	Phenol, 2,4,6-tris(1-methylethyl)-	88.63	88.51	89.75				5	Butylated hydroxy toluene	87.63	87.09	92.45				6	Phenol, 2,4-bis(1,1-dimethylethyl)-6-...	86.72	86.64	87.44				7	(5S,6S)-6-(cyclohepta-2,4,6-trien-1-yl)-...	86.26	85.78	90.55				8	4-Hydroxy-3,5-diisopropylbenzaldehy...	81.40	81.34	81.93				9	Ethanone, 1-(5,6,7,8-tetrahydro-2,8,8-...	78.30	78.25	78.77		
RT [min]	#	Match	Score	HQI	R.H...	Com...	N																																																																																																																																																																			
3.8140	1	Trimethyl-phosphate	89.17	89.12	89.61																																																																																																																																																																					
3.8952	1	alpha-Pinene	91.44	91.41	91.73																																																																																																																																																																					
6.8448	1	Biphenyl	97.09	96.96	98.27																																																																																																																																																																					
6.9482	1	C-Phenyl-C-formyl-N-(2-furyl)nitron	90.20	89.52	96.24																																																																																																																																																																					
6.9492	1	Diphenylether	90.85	90.77	91.50																																																																																																																																																																					
7.0295	1	Diphenylether	76.30	75.40	84.41																																																																																																																																																																					
7.0370	1	1,1'-Biphenyl, 3-ethoxy-	68.56	68.19	71.88																																																																																																																																																																					
7.1110	1	Biphenyl	55.72	51.79	91.15																																																																																																																																																																					
7.2828	1	Phenol, 3-(1,1-dimethylethyl)-4-metho...	93.02	92.97	93.49																																																																																																																																																																					
7.3261	1	3-tert-Butyl-4-hydroxyanisol	95.83	95.80	96.11																																																																																																																																																																					
7.4200	1	Dimethylisophthalate	90.00	89.93	90.65																																																																																																																																																																					
7.4282	1	Butylated hydroxytoluene	93.35	93.25	94.26																																																																																																																																																																					
	2	Phenol, 2,6-bis(1,1-dimethylethyl)-4-...	90.14	90.05	90.96																																																																																																																																																																					
	3	Phenol, 2,6-bis(1,1-dimethylethyl)-4-...	90.09	89.98	91.07																																																																																																																																																																					
	4	Phenol, 2,4,6-tris(1-methylethyl)-	88.63	88.51	89.75																																																																																																																																																																					
	5	Butylated hydroxy toluene	87.63	87.09	92.45																																																																																																																																																																					
	6	Phenol, 2,4-bis(1,1-dimethylethyl)-6-...	86.72	86.64	87.44																																																																																																																																																																					
	7	(5S,6S)-6-(cyclohepta-2,4,6-trien-1-yl)-...	86.26	85.78	90.55																																																																																																																																																																					
	8	4-Hydroxy-3,5-diisopropylbenzaldehy...	81.40	81.34	81.93																																																																																																																																																																					
	9	Ethanone, 1-(5,6,7,8-tetrahydro-2,8,8-...	78.30	78.25	78.77																																																																																																																																																																					

	Action	Result
5	<p>To export results to ReportIt, use the Transfer to: toolbar. Click on ReportIt.</p> <p>A window will pop-up. Select the a report template, GC_Expert_Portrait.</p> <p>Click OK.</p>	

	Action	Result
6		<div data-bbox="825 321 1839 760"> </div> <p>Note: If this is the first time using the GC Expert reporting feature, templates will need to be added ahead of this step. The following steps should be performed using GC Expert.</p> <ul style="list-style-type: none"> • File > Edit Report Templates • Choose Add button in the pop-up window. • Select the files, GC_Expert Landscape and GC_Expert Portrait (located "C:\Users\Public\Documents\Wiley\KnowItAll\Report Templates\GC Expert") • Click Open. • Click Close in the Report Templates pop-up window.

Action		Result								
7	Go back to GC Expert .									
	Expand a few components' reference hits, and add some comments									
	Right-mouse click the Components table and select Copy All Component Information . One can then paste this information to Excel or Word, etc.									
		1	RT [min]	#	Match	Score	HQI	R.HQI	Notes	Area %
		2	1.3211		1 (E)-1,2-bis	81.69	81.49	83.55		0.01
		3	1.357		1 Trichlorom	76.85	75.94	85.11		4.05
		4	1.4468		1 (E)-1,2-bis	94.86	94.52	97.85		0.05
		5	1.9667		1 Cyclopent	90.92	90.8	92.06		0.07
		6	3.814		1 Trimethyl p	88.17	87.8	91.46		4.32
		7	3.8972		1 alpha-Pine	91.35	91.28	91.99		3.73
		8	6.8449		1 Biphenyl	96.72	96.56	98.16		5.59
		9	6.9482		1 C-Phenyl-C	90.27	89.61	96.29		8.26
		10	6.9492		1 Diphenyltet	91.11	91.04	91.75		8.32
		11	7.0295		1 Diphenyltet	87.03	86.92	88.03	This is righ	0.04
		12			2 DIPHENYL	85.7	84.87	93.22		
		13			3 C-Phenyl-C	84.56	83.32	95.75		
		14			4 2-Phenyl-p	84.5	84.18	87.39		
		15			5 Biphenyltol	82.6	82.52	83.33		
		16			6 8-Methyl-1	81.62	80.31	93.42		
		17			7 5-Methyl-1	81.57	80.26	93.36		
		18			8 1,2-Dihydr	80.56	79.27	92.21		
		19			9 1H,3H-Naj	80.43	80.31	81.46		
		20			10 1-Naphtha	80.15	78.86	91.74		
		21	7.1114		1 Biphenyl	72.39	70.33	90.86		0
		22	7.2828		1 Phenol, 3-i	93.06	93.02	93.43		0.06
		23	7.3261		1 3-tert-Buty	95.51	95.47	95.86		3.88
		24	7.4202		1 Dimethylis	90.19	90.14	90.67		13.89
		25	7.4282		1 Butylated l	93.18	93.08	94.04		6.35
		26	7.5162		1 Dimethylis	76.71	75.93	83.71		0.01
		27	7.8488		1 Diethylpht	94.56	94.43	95.78		9.19
		28	9.083		1 Caffeine	90.93	90.78	92.26		3.34
		29	9.3082		1 Hexadecar	92.73	92.51	94.73		3.1
		30	9.4719		1 1,4-Dibuty	95.38	95.27	96.35		6.57
		31	9.473		1 Dibutyl pht	93.81	93.76	94.24		7.52
		32	10.0283		1 Drometriz	95.72	95.7	95.94		0.97
	33	10.3739		1 Phenol, 4,	91.76	91.63	92.96		0.85	
	34	10.4934		1 m-Terpher	97.01	96.93	97.72		9.83	

How to Perform GC-MS Analysis by Peak Picking

Purpose

These exercises demonstrate how to use KnowItAll GC Expert to analyze GC-MS data by peak picking.

Objectives

These exercises will teach you:

- How to use KnowItAll GC Expert to analyze GC-MS using peak picking.
 - How to generate reports.
-

Background

GC-MS data can be analyzed by peak picking without going through deconvolution. Novel compounds can be identified and structural characteristics can be deduced from applying the MS Adaptive search that uses fragmentation and structural data to propose likely structural details of the unknown.

Training Files Used in This Lesson


- C:\Users\Public\Documents\Wiley\KnowItAll\Samples\GC-MS\Barbiturate GC-MS.d



KnowItAll Applications Used

- KnowItAll GC Expert

KnowItAll GC Expert defaults to perform deconvolution. However, one can also use the **Peak Picking** option to analyze GC-MS data, especially if the GC separates components well enough.

Example: Peak Picking

	Action	Result
1	<p>Click Open Raw GC-MS Data File button.</p> <p>Navigate to “C:\Users\Public\Documents\Wiley\KnowItAll\Samples\GC-MS”</p> <p>Select the file: Barbiturate GC-MS.d</p> <p>And click Open.</p>	 <p>The screenshot shows the 'Open' file dialog box with the 'Look in' field set to 'GC-MS'. The file list contains several files, with 'Barbiturate GC-MS.d' selected. Below the file list, the 'File name' field is set to 'Barbiturate GC-MS' and the 'Files of type' is set to 'All Files (*.*)'. The 'Open' button is highlighted. Below the file list, there is a section for 'Imported spectrum is' and a chromatogram plot. The plot shows a baseline with several sharp peaks, labeled 'BARBITURATE MIX'. The x-axis is labeled 'min' and ranges from 1 to 6. The y-axis is labeled 'Data is 2D, only the TIC is shown.'</p>

	Action	Result
2	Switch Analysis Method to Peak Picking .	
3	<p>Reduce the Sensitivity to Low so that only main peaks are picked.</p> <p>Select the peak at 3.8917 min.</p> <p>Hold the mouse cursor over the solid square denoting the selected peak in the GC panel.</p> <ul style="list-style-type: none"> In the Raw Spectrum panel, B displays the MS spectrum of the selected peak C, in both the GC and Database Match panels, displays the selected peaks <p>In the Database Match panel, D indicates the reference spectrum match to the selected spectrum indicated by C.</p>	

Manually Analyze GC-MS Data

How to manually analyze GC-MS data

Purpose

These exercises demonstrate how to use KnowItAll GC Expert to manually analyze GC-MS data.

Objectives

These exercises will teach you:

- How to use KnowItAll GC Expert to manually analyze GC-MS data.
-

Background

KnowItAll GC Expert allows one to examine GC-MS data and perform spectral subtraction. Matches to the subtracted spectrum are searched for against the reference data.


Training Files Used in This Lesson



- C:\Users\Public\Documents\Wiley\KnowItAll\Samples\GC-MS\Barbiturate GC-MS.d

KnowItAll Applications Used

- KnowItAll GC Expert

Example: Manual Analysis

	Action	Result
1	<p>In GC Expert, click on the Open Raw GC-MS Data File button.</p> <p>Navigate to: "C:\Users\Public\Documents\Wiley\KnowItAll\Samples\GC-MS"</p> <p>Select the file: Barbiturate GC-MS.d</p> <p>Click Open.</p>	

	Action	Result																																																																													
2	Switch Analysis Method to Manual .																																																																														
3	<p>Use the Average Range Bar (A) to define a selected region of the MS chromatogram.</p> <ul style="list-style-type: none">Click and drag the mouse cursor along the Average Range Bar.You can select as many ranges as required. <p>Use the Background Range Bar (B) to define the selected background region of the chromatogram.</p> <ul style="list-style-type: none">Click and drag the mouse cursor along the Background Range Bar. <p>You can select as many ranges as required.</p>	<p>The MS spectrum with the subtracted background (C) and the corresponding matched reference MS spectrum (D) is displayed in the bottom panel.</p>  <table><thead><tr><th>RT [min]</th><th># Match</th><th>Score</th><th>HQI</th><th>R.H.L.</th><th>Notes</th><th>Are...</th></tr></thead><tbody><tr><td>2.2948...</td><td>1</td><td>94.64</td><td>94.64</td><td>94.70</td><td></td><td></td></tr><tr><td>2</td><td>Barbitone</td><td>88.47</td><td>87.91</td><td>93.47</td><td></td><td></td></tr><tr><td>3</td><td>2,4,6-(TH,3H,5H)-Pyrimidinotri...</td><td>88.00</td><td>86.72</td><td>99.46</td><td></td><td></td></tr><tr><td>4</td><td>Barbitone</td><td>86.43</td><td>86.10</td><td>89.44</td><td></td><td></td></tr><tr><td>5</td><td>Hawethal</td><td>85.82</td><td>85.29</td><td>90.60</td><td></td><td></td></tr><tr><td>6</td><td>2,4,6-(TH,3H,5H)-Pyrimidinotri...</td><td>85.44</td><td>85.05</td><td>88.92</td><td></td><td></td></tr><tr><td>7</td><td>2,4,6-(TH,3H,5H)-Pyrimidinotri...</td><td>85.23</td><td>84.90</td><td>88.25</td><td></td><td></td></tr><tr><td>8</td><td>Butobarbital</td><td>84.96</td><td>84.96</td><td>84.96</td><td></td><td></td></tr><tr><td>9</td><td>Pentobarbitalone</td><td>81.64</td><td>81.33</td><td>84.49</td><td></td><td></td></tr><tr><td>10</td><td>Probarbital</td><td>81.39</td><td>81.38</td><td>81.44</td><td></td><td></td></tr></tbody></table>	RT [min]	# Match	Score	HQI	R.H.L.	Notes	Are...	2.2948...	1	94.64	94.64	94.70			2	Barbitone	88.47	87.91	93.47			3	2,4,6-(TH,3H,5H)-Pyrimidinotri...	88.00	86.72	99.46			4	Barbitone	86.43	86.10	89.44			5	Hawethal	85.82	85.29	90.60			6	2,4,6-(TH,3H,5H)-Pyrimidinotri...	85.44	85.05	88.92			7	2,4,6-(TH,3H,5H)-Pyrimidinotri...	85.23	84.90	88.25			8	Butobarbital	84.96	84.96	84.96			9	Pentobarbitalone	81.64	81.33	84.49			10	Probarbital	81.39	81.38	81.44		
RT [min]	# Match	Score	HQI	R.H.L.	Notes	Are...																																																																									
2.2948...	1	94.64	94.64	94.70																																																																											
2	Barbitone	88.47	87.91	93.47																																																																											
3	2,4,6-(TH,3H,5H)-Pyrimidinotri...	88.00	86.72	99.46																																																																											
4	Barbitone	86.43	86.10	89.44																																																																											
5	Hawethal	85.82	85.29	90.60																																																																											
6	2,4,6-(TH,3H,5H)-Pyrimidinotri...	85.44	85.05	88.92																																																																											
7	2,4,6-(TH,3H,5H)-Pyrimidinotri...	85.23	84.90	88.25																																																																											
8	Butobarbital	84.96	84.96	84.96																																																																											
9	Pentobarbitalone	81.64	81.33	84.49																																																																											
10	Probarbital	81.39	81.38	81.44																																																																											