

KnowItAll[®] Informatics Training

NMR Processing Tools

Process NMR

How to Import and Process 1D NMR Files

Purpose

This exercise demonstrates how to use the KnowItAll Informatics System's ProcessIt NMR application to import and process 1D NMR files.

Objectives

This exercise will teach you:

- How to use ProcessIt NMR tools to improve spectral appearance and correct experimental artifacts
 - How to create and use processing macros
-

Background

You can use the ProcessIt application to import raw data files from major NMR instrument vendors and processed data formats. The application can then be used to process the files and is useful for improving spectral appearance and correcting experimental artifacts.



Training Files Used in This Lesson

- C:\Users\Public\Public Documents\Wiley\KnowItAll\Samples\NMR\Bruker TopSpin\C14H10ClNO3\C14H10ClNO3 H1\fid




KnowItAll Applications Used



- ProcessIt


Open a raw NMR file in ProcessIt



	Action	Result
1	<p>Open the ProcessIt application by clicking its icon, typically found in the Spectral Processing toolbox.</p>  ProcessIt	
2	<p>In ProcessIt, click Open Data File button. Navigate to "C:\Users\Public\Public Documents\Wiley\KnowItAll\Samples\NMR\Bruker TopSpin\C14H10ClNO3\C14H10ClNO3 H1\fid".</p> <p>Note: ProcessIt application can import many different NMR formats and nuclei, and can open either the raw (e.g., fid) or processed file.</p>	<p>The file is displayed in the main ProcessIt window:</p> 



Methods to Transform the File


	Action	Result
1	<p>Fourier Transform: If you wish to apply a Fourier Transform (FT) without additional preprocessing, begin here (<i>otherwise move to step 4</i>).</p> <p>Choose Process > Transforms > Fourier Transform, or click the Fourier Transform icon on the Process Toolbar.</p> 	<p>The Fourier Transform dialog appears with a preview of the transformed proton spectrum behind it:</p> 
2	<p>Click Replace on the Fourier Transform dialog to accept the previewed transformation using the Default algorithm.</p> <p>Note: Four different algorithms are available in the Fourier Transform dialog: Default, Complex, Real and Bruker.</p>	<p>The transformed proton spectrum is displayed:</p> 

	Action	Result
3	<p>Note: The next steps describe how to apply pre-processing to a raw file before performing the transformation. If you are ready to process your transformed file, move to the next section "Methods to Process the Spectrum".</p>	
4	<p>Begin by reopening the file in Step 1. Choose File > Open Data File or click the Open Data File icon on the Standard Toolbar.</p> 	<p>The file is displayed in the main ProcessIt window:</p> 



5 Zero Filling: Choose Process > Zero Fill .	<p>The Zero Filling dialog box opens:</p>  <p>The software proposes a new number of points for the spectrum based on the original (experimental) count.</p>
--	--



	Action	Result
6	Click Replace .	<p>The software executes the Zero Filling command:</p> 
7	<p>Window Functions: Choose Process > Window Functions.</p>	<p>The Window Functions dialog box opens:</p>  <p>Window Functions are used extensively to increase resolution and decrease noise in the spectrum.</p>



	Action	Result
8	<p>Select Exponential under Algorithm. Check Perform Auto DC to apply a DC Offset correction to the spectrum. Select Default under FT Preview options to apply a Fourier Transform to the processed spectrum.</p> <p>Note: The line broadening parameter (LB) can also be modified if desired.</p>	
	<p>Click Preview to view the resulting transformed spectrum.</p>	<p>A preview of the transformed spectrum according to the options in the Window Functions dialog is displayed:</p> 



	Action	Result
9	Click Replace to accept the previewed Fourier Transform spectrum.	<p>The transformed spectrum is displayed in ProcessIt:</p>  <p>The screenshot displays the ProcessIt software window. The title bar reads 'ProcessIt'. The main plot area shows a 1H NMR spectrum for the sample 'C14H10ClNO3 H1'. The x-axis is labeled 'ppm' and ranges from 15 to -1. The y-axis represents intensity, ranging from 0 to 6e+08. A prominent peak is visible at approximately 7.2 ppm. Other smaller peaks are present in the aromatic region (6.5-8.5 ppm) and the aliphatic region (1.5-4.5 ppm). The software interface includes a toolbar with various processing tools and a status bar at the bottom showing the sample name 'C14H10ClNO3 H1'.</p>



Methods to process the spectrum


	Action	Result
1	<p>Phasing: Choose Process > Phasing.</p> <p>Note: Manual phasing can be performed by clicking and dragging with either the left (PH0) or right (PH1) mouse button.</p>	<p>The Phase Correction dialog box opens with a preview of the phased spectrum applying the phase angles in Current phase angles (degrees).</p>  <p>Three different automatic phasing methods are available in the Phase Correction dialog: GoodLook, PAMPAS (Phase Angle Measurement from Peak Area), and APSL (Automatic 1D Phase Correction).</p>
2	<p>Retain GoodLook as the automatic phasing method and click Replace.</p>	<p>The phasing for the spectrum is improved upon applying the Phase Correction function:</p> 



	Action	Result
3	<p>Baseline Correction: Choose Process > Baseline Correction.</p>	<p>The Baseline Correction dialog box opens and a preview of the processed spectrum is displayed:</p>  <p>Spline, linear, polynomial and FaceLift algorithms can be applied to baseline points selected either automatically or manually. As a technique is selected, the resulting spectrum is previewed.</p>
4	<p>Choose Linear as the baseline technique. Click Replace.</p>	<p>The Baseline Correction function is applied:</p> 



	Action	Result
5	Smoothing: Choose Analysis > Smooth Spectrum .	The Smooth NMR Spectrum popup window opens: 
6	Enter 0.0003 as value in the Filter width (ppm) box, then click Replace .	The Smooth NMR Spectrum popup closes and the smoothed spectrum is displayed: 



	Action	Result
7	<p>Peak Picking: Choose Process > Peak Picking.</p>	<p>The Peak picking dialogue appears. On the spectrum, a preview for the peaks are displayed as vertical lines using Automatic Peak Picking method:</p> 
8	<p>Enter 4 as the value for Minimum Intensity (%) and select Preview to preview the result of the peak picking. Choose Replace to accept the changes.</p>	<p>The peaks appear on the spectrum:</p> 



	Action	Result																								
9	<p>On the Peak picking dialog, select the Manual tab.</p>	<p>The Manual Peak Picking table is visible and the Peak Bar is displayed on the spectrum:</p>  <table border="1" data-bbox="856 500 1207 706"> <thead> <tr> <th>Peak [ppm]</th> <th>Intensit</th> </tr> </thead> <tbody> <tr><td>-7.4652e-05</td><td>29781509.11</td></tr> <tr><td>1.50</td><td>27231358.88</td></tr> <tr><td>1.50</td><td>30922249.84</td></tr> <tr><td>1.50</td><td>31754208.01</td></tr> <tr><td>1.51</td><td>32616367.34</td></tr> <tr><td>1.51</td><td>32867922.33</td></tr> <tr><td>1.51</td><td>32948125.61</td></tr> <tr><td>1.51</td><td>32938471.09</td></tr> <tr><td>1.51</td><td>33043031.56</td></tr> <tr><td>1.51</td><td>33030814.27</td></tr> <tr><td>1.51</td><td>32742630.35</td></tr> </tbody> </table>	Peak [ppm]	Intensit	-7.4652e-05	29781509.11	1.50	27231358.88	1.50	30922249.84	1.50	31754208.01	1.51	32616367.34	1.51	32867922.33	1.51	32948125.61	1.51	32938471.09	1.51	33043031.56	1.51	33030814.27	1.51	32742630.35
Peak [ppm]	Intensit																									
-7.4652e-05	29781509.11																									
1.50	27231358.88																									
1.50	30922249.84																									
1.50	31754208.01																									
1.51	32616367.34																									
1.51	32867922.33																									
1.51	32948125.61																									
1.51	32938471.09																									
1.51	33043031.56																									
1.51	33030814.27																									
1.51	32742630.35																									
10	<p>Click on the first row in the Manual Peak Picking table. Press CTRL on the keyboard and click each row with ppm value less than 3.91 ppm so that they become highlighted. Press Delete button on keyboard, followed by clicking Replace on the dialog.</p> <p>Note: Any row can be deleted one at a time by pressing Delete button on keyboard.</p>	<p>The peaks less than 3.91 ppm are removed from the Manual Peak Picking table and spectrum:</p>  <table border="1" data-bbox="856 1052 1207 1242"> <thead> <tr> <th>Peak [ppm]</th> <th>Intensit</th> </tr> </thead> <tbody> <tr><td>3.91</td><td>39430165.53</td></tr> <tr><td>3.92</td><td>55323671.48</td></tr> <tr><td>3.93</td><td>43595044.52</td></tr> <tr><td>4.38</td><td>52446780.06</td></tr> <tr><td>4.39</td><td>66238649.92</td></tr> <tr><td>4.41</td><td>49905938.73</td></tr> <tr><td>7.19</td><td>664559214.93</td></tr> <tr><td>7.77</td><td>58257232.23</td></tr> <tr><td>7.79</td><td>91296597.38</td></tr> <tr><td>7.81</td><td>33155959.49</td></tr> <tr><td>7.81</td><td>33008482.23</td></tr> </tbody> </table>	Peak [ppm]	Intensit	3.91	39430165.53	3.92	55323671.48	3.93	43595044.52	4.38	52446780.06	4.39	66238649.92	4.41	49905938.73	7.19	664559214.93	7.77	58257232.23	7.79	91296597.38	7.81	33155959.49	7.81	33008482.23
Peak [ppm]	Intensit																									
3.91	39430165.53																									
3.92	55323671.48																									
3.93	43595044.52																									
4.38	52446780.06																									
4.39	66238649.92																									
4.41	49905938.73																									
7.19	664559214.93																									
7.77	58257232.23																									
7.79	91296597.38																									
7.81	33155959.49																									
7.81	33008482.23																									




	Action	Result																								
11	On the spectrum, click on the Peak Bar at ~1.50 ppm. Then click Replace on the Peak picking dialog.	A new peak is added to the Manual Peak Picking table. A vertical line representing the new peak is displayed on the spectrum with a peak label of the corresponding ppm:  <table border="1" data-bbox="850 527 1165 755"><thead><tr><th>Peak (ppm)</th><th>Intensity</th></tr></thead><tbody><tr><td>1.51</td><td>32376268.97</td></tr><tr><td>3.91</td><td>39430165.53</td></tr><tr><td>3.92</td><td>55323671.48</td></tr><tr><td>3.93</td><td>43556044.52</td></tr><tr><td>4.38</td><td>52446780.06</td></tr><tr><td>4.39</td><td>66238649.92</td></tr><tr><td>4.41</td><td>49905938.73</td></tr><tr><td>7.19</td><td>664559214.93</td></tr><tr><td>7.77</td><td>58257232.23</td></tr><tr><td>7.79</td><td>91296597.38</td></tr><tr><td>7.81</td><td>33155959.49</td></tr></tbody></table>	Peak (ppm)	Intensity	1.51	32376268.97	3.91	39430165.53	3.92	55323671.48	3.93	43556044.52	4.38	52446780.06	4.39	66238649.92	4.41	49905938.73	7.19	664559214.93	7.77	58257232.23	7.79	91296597.38	7.81	33155959.49
Peak (ppm)	Intensity																									
1.51	32376268.97																									
3.91	39430165.53																									
3.92	55323671.48																									
3.93	43556044.52																									
4.38	52446780.06																									
4.39	66238649.92																									
4.41	49905938.73																									
7.19	664559214.93																									
7.77	58257232.23																									
7.79	91296597.38																									
7.81	33155959.49																									
12	Choose Close to close the Peak picking dialog.	The spectrum is displayed with peak labels for the picked peaks: 																								





	Action	Result
13	<p>Reference Spectrum: Use this tool to reference the spectrum to a solvent or internal standard.</p> <p>Choose Analysis > Set Reference.</p>	<p>The Set Reference popup window appears:</p> 
14	<p>On the spectrum, double click on the peak close to 7.2 ppm to set it as the reference point. Change New Value (ppm) to 7.26. Click Apply.</p>	<p>Upon clicking Apply, the Current Value (ppm) becomes 7.26. On the spectrum, all of the peaks move according to the reference value. The spectrum reference point is shifted to 7.26 ppm:</p> 






	Action	Result
15	Click OK to close the dialog.	<p>The Set Reference dialog closes:</p> 
16	<p>Automatic Integration: Choose Analysis > Peak Area/Integration.</p> <p>Note: Two methods for automatic integration are available: Perpendicular Drop and Tangential Skim. By default, Tangential Skim is preselected.</p>	<p>The Peak Area/Integration Table popup appears:</p> 











	Action	Result																																																								
17	<p>Choose the radio button for Perpendicular Drop, then click Automatic Integration.</p> <p>The Automatic Integration dialog box allows you to specify Noise Threshold and Minimum Distance between Peaks, and to Perform Baseline Correction before determining the integration.</p> <p>Note: Perpendicular Drop is suggested for NMR automatic integration.</p>	<p>The Automatic Integration dialog opens:</p>  <p>If baseline correction is applied, an automatic baseline correction will be performed prior to peak integration. If it is unchecked, the raw data will be integrated.</p>																																																								
18	<p>Click Preview then Replace.</p> <p>Note: Any cell in the Peak Area/Integration Table can be modified by double clicking directly on the cell, or by selecting a row and clicking Edit Range, Edit Raw Peak Area, or Edit Integration Value.</p> <p>Note: Any row can be deleted by clicking Delete Row, or all rows can be deleted by clicking Clear All.</p>	<p>The Peak Area/Integration Table becomes populated with integral information: left x-value (X₁), right x-value (X₂), Raw Peak Area, Integration, Peak Height, and full width at half maximum (FWHM). The integrals are visible on the spectrum:</p>  <table border="1" data-bbox="695 1013 1331 1170"> <thead> <tr> <th>ID</th> <th>X₁ [ppm]</th> <th>X₂ [ppm]</th> <th>Raw Peak Area</th> <th>Integration</th> <th>Peak Height</th> <th>FWHM [ppm]</th> </tr> </thead> <tbody> <tr> <td>1</td> <td>8.696</td> <td>8.669</td> <td>296188.1916</td> <td>2.6485</td> <td>35012504.391</td> <td>0.018</td> </tr> <tr> <td>2</td> <td>8.628</td> <td>8.621</td> <td>111833.7938</td> <td>1.0000</td> <td>36032465.336</td> <td>0.003</td> </tr> <tr> <td>3</td> <td>8.541</td> <td>8.511</td> <td>334952.7666</td> <td>2.9951</td> <td>62998785.309</td> <td>0.020</td> </tr> <tr> <td>4</td> <td>7.863</td> <td>7.830</td> <td>412971.1325</td> <td>3.6927</td> <td>91296597.385</td> <td>0.021</td> </tr> <tr> <td>5</td> <td>7.272</td> <td>7.250</td> <td>1371003.4421</td> <td>12.2593</td> <td>664559214.925</td> <td>0.001</td> </tr> <tr> <td>6</td> <td>4.480</td> <td>4.439</td> <td>695635.1622</td> <td>6.2203</td> <td>66238649.921</td> <td>0.016</td> </tr> <tr> <td>7</td> <td>4.009</td> <td>3.967</td> <td>700012.8139</td> <td>6.2594</td> <td>55323671.483</td> <td>0.029</td> </tr> </tbody> </table>	ID	X ₁ [ppm]	X ₂ [ppm]	Raw Peak Area	Integration	Peak Height	FWHM [ppm]	1	8.696	8.669	296188.1916	2.6485	35012504.391	0.018	2	8.628	8.621	111833.7938	1.0000	36032465.336	0.003	3	8.541	8.511	334952.7666	2.9951	62998785.309	0.020	4	7.863	7.830	412971.1325	3.6927	91296597.385	0.021	5	7.272	7.250	1371003.4421	12.2593	664559214.925	0.001	6	4.480	4.439	695635.1622	6.2203	66238649.921	0.016	7	4.009	3.967	700012.8139	6.2594	55323671.483	0.029
ID	X ₁ [ppm]	X ₂ [ppm]	Raw Peak Area	Integration	Peak Height	FWHM [ppm]																																																				
1	8.696	8.669	296188.1916	2.6485	35012504.391	0.018																																																				
2	8.628	8.621	111833.7938	1.0000	36032465.336	0.003																																																				
3	8.541	8.511	334952.7666	2.9951	62998785.309	0.020																																																				
4	7.863	7.830	412971.1325	3.6927	91296597.385	0.021																																																				
5	7.272	7.250	1371003.4421	12.2593	664559214.925	0.001																																																				
6	4.480	4.439	695635.1622	6.2203	66238649.921	0.016																																																				
7	4.009	3.967	700012.8139	6.2594	55323671.483	0.029																																																				

	Action	Result																												
19	<p>Note: Integral ranges can be preloaded into the Peak Area/Integration Table dialog.</p> <p>To use this feature, click Load Ranges and import a CSV file that has the format X1, X2. For example:</p> <p>X1,X2 8.0,7.6 4.6,4.2 4.1,3.7</p>	<p>The integrals are calculated for the imported regions:</p>  <table border="1" data-bbox="695 505 1289 581"> <thead> <tr> <th>ID</th> <th>X₁ [ppm]</th> <th>X₂ [ppm]</th> <th>Raw Peak Area</th> <th>Integration</th> <th>Peak Height</th> <th>FWHM [ppm]</th> </tr> </thead> <tbody> <tr> <td>1</td> <td>8</td> <td>7.600</td> <td>751929.0629</td> <td>1.0088</td> <td>91296597.385</td> <td>0.021</td> </tr> <tr> <td>2</td> <td>4.600</td> <td>4.200</td> <td>751891.6205</td> <td>1.0088</td> <td>66238649.921</td> <td>0.016</td> </tr> <tr> <td>3</td> <td>4.100</td> <td>3.700</td> <td>745333.8324</td> <td>1.0000</td> <td>55323671.483</td> <td>0.029</td> </tr> </tbody> </table>	ID	X ₁ [ppm]	X ₂ [ppm]	Raw Peak Area	Integration	Peak Height	FWHM [ppm]	1	8	7.600	751929.0629	1.0088	91296597.385	0.021	2	4.600	4.200	751891.6205	1.0088	66238649.921	0.016	3	4.100	3.700	745333.8324	1.0000	55323671.483	0.029
ID	X ₁ [ppm]	X ₂ [ppm]	Raw Peak Area	Integration	Peak Height	FWHM [ppm]																								
1	8	7.600	751929.0629	1.0088	91296597.385	0.021																								
2	4.600	4.200	751891.6205	1.0088	66238649.921	0.016																								
3	4.100	3.700	745333.8324	1.0000	55323671.483	0.029																								
20	<p>Click Close to exit the Peak Area/Integration Table dialog.</p>	<p>The Peak Area/Integration Table dialog is closed:</p> 																												




	Action	Result
21	<p>Manual Integration: Choose Analysis > Peak Area/Integration. To clear any previous integrals, click Clear All.</p>	<p>The Peak Area/Integration Table popup launches. The integrals are cleared from the table and spectrum:</p> 
22	<p>Move or resize the Peak/Area Integration Table popup so that the peaks on the spectrum become visible. Right click on the spectrum and select Horizontal Zoom Mode. Click and drag the mouse cursor from ~ 7 to 9 ppm. Let go of the mouse cursor.</p>	<p>The mouse cursor changes to a zoom tool () and the spectrum is focused on the region from ~ 7 to 9 ppm.</p> 



	Action	Result
23	Select the Peak Area/Integration Table icon located on the Process toolbar to reactivate the integration tool. 	Clicking the Peak Area/Integration Table icon changes the mouse cursor from the zoom tool () to the integration tool ():
24	Click and drag the mouse button across a multiplet. Release the mouse button at the end of the multiplet. Repeat for each multiplet in visible in the displayed region. The multiplets are: (1) 8.67 – 8.70 ppm, (2) 8.62 – 8.65 ppm, (3) 8.51 – 8.54, (4) 7.83 – 7.91 ppm. Note: Dragging functionality is bidirectional (left to right, or right to left). Note: Do not integrate the peak at 7.26 ppm which represents the ¹ H-NMR solvent signal.	The integration value for each multiplet is added as a new row in the Peak Area/Integration Table . The integral also display on the spectrum.  <p>By default, the integral values are normalized to 1.0 for the smallest integral. To change the normalization to a different value (e.g., to 2.0), select a row in the table and click Edit Integration Value. Change the value in the proceeding Integration Value popup and click OK. The other integrals will adjust accordingly.</p>



	Action	Result																																			
25	<p>Right click on the spectrum and click pan mode. Click and hold the left mouse button. Drag the spectrum to the left to move the spectrum to a lower ppm region. Release the mouse button when more multiplets are visible (~3.2 ppm to 5 ppm).</p>	<p>Upon selecting pan mode, the corresponding mouse cursor is displayed () which allows the spectrum to move to the right or left.</p>  <table border="1" data-bbox="808 581 1333 665"> <thead> <tr> <th>ID</th> <th>X₁ [ppm]</th> <th>X₂ [ppm]</th> <th>Raw Peak Area</th> <th>Integration</th> <th>Peak Height</th> <th>FWHM [ppm]</th> </tr> </thead> <tbody> <tr> <td>1</td> <td>7.900</td> <td>7.830</td> <td>619865.0622</td> <td>2.2525</td> <td>90255306.341</td> <td>0.002</td> </tr> <tr> <td>2</td> <td>8.545</td> <td>8.509</td> <td>314796.6795</td> <td>1.1461</td> <td>62414455.383</td> <td>0.020</td> </tr> <tr> <td>3</td> <td>8.658</td> <td>8.611</td> <td>337745.9263</td> <td>1.2297</td> <td>35627850.943</td> <td>0.021</td> </tr> <tr> <td>4</td> <td>8.702</td> <td>8.672</td> <td>274661.2948</td> <td>1.0000</td> <td>32371628.412</td> <td>0.018</td> </tr> </tbody> </table>	ID	X ₁ [ppm]	X ₂ [ppm]	Raw Peak Area	Integration	Peak Height	FWHM [ppm]	1	7.900	7.830	619865.0622	2.2525	90255306.341	0.002	2	8.545	8.509	314796.6795	1.1461	62414455.383	0.020	3	8.658	8.611	337745.9263	1.2297	35627850.943	0.021	4	8.702	8.672	274661.2948	1.0000	32371628.412	0.018
ID	X ₁ [ppm]	X ₂ [ppm]	Raw Peak Area	Integration	Peak Height	FWHM [ppm]																															
1	7.900	7.830	619865.0622	2.2525	90255306.341	0.002																															
2	8.545	8.509	314796.6795	1.1461	62414455.383	0.020																															
3	8.658	8.611	337745.9263	1.2297	35627850.943	0.021																															
4	8.702	8.672	274661.2948	1.0000	32371628.412	0.018																															
26	<p>To deactivate pan mode and activate the integration tool, reselect the Peak Area/Integration Table icon located on the Process toolbar. </p> <p>Integrate the two multiplets in this region by dragging the mouse button across each multiplet.</p>	<p>Selecting the Peak Area/Integration Table icon changes the mouse cursor back to the integration tool () ,changing the functionality of the mouse cursor. Two new integrals are added to the spectrum:</p>  <table border="1" data-bbox="808 1128 1333 1213"> <thead> <tr> <th>ID</th> <th>X₁ [ppm]</th> <th>X₂ [ppm]</th> <th>Raw Peak Area</th> <th>Integration</th> <th>Peak Height</th> <th>FWHM [ppm]</th> </tr> </thead> <tbody> <tr> <td>3</td> <td>8.658</td> <td>8.611</td> <td>337745.9263</td> <td>1.2297</td> <td>35627850.943</td> <td>0.021</td> </tr> <tr> <td>4</td> <td>8.702</td> <td>8.672</td> <td>274661.2948</td> <td>1.0000</td> <td>32371628.412</td> <td>0.018</td> </tr> <tr> <td>5</td> <td>4.483</td> <td>4.437</td> <td>643625.7707</td> <td>2.3433</td> <td>64896811.967</td> <td>0.027</td> </tr> <tr> <td>6</td> <td>4.008</td> <td>3.962</td> <td>619419.5592</td> <td>2.2552</td> <td>53389709.500</td> <td>0.024</td> </tr> </tbody> </table>	ID	X ₁ [ppm]	X ₂ [ppm]	Raw Peak Area	Integration	Peak Height	FWHM [ppm]	3	8.658	8.611	337745.9263	1.2297	35627850.943	0.021	4	8.702	8.672	274661.2948	1.0000	32371628.412	0.018	5	4.483	4.437	643625.7707	2.3433	64896811.967	0.027	6	4.008	3.962	619419.5592	2.2552	53389709.500	0.024
ID	X ₁ [ppm]	X ₂ [ppm]	Raw Peak Area	Integration	Peak Height	FWHM [ppm]																															
3	8.658	8.611	337745.9263	1.2297	35627850.943	0.021																															
4	8.702	8.672	274661.2948	1.0000	32371628.412	0.018																															
5	4.483	4.437	643625.7707	2.3433	64896811.967	0.027																															
6	4.008	3.962	619419.5592	2.2552	53389709.500	0.024																															


	Action	Result
27	<p>Note: On the spectrum, any integral can be resized. To resize an integral, click on the edge of the integral region and drag its edge to left or right along the x-axis.</p>	<p>When the mouse cursor is positioned on the edge of the integral region, the cursor changes from the integral cursor () to display an arrow pointing right () or left (). The example integral below becomes wider upon dragging the integral edge to the right:</p> <p>Before dragging right integral edge:  After dragging integral edge to the right: </p>
28	<p>Note: On the spectrum, any integral can be shifted. To relocate an integral, click inside of the integral region and drag it to the left or right on the x-axis.</p>	<p>When the mouse cursor is positioned inside the integral region, the cursor changes from the integral cursor () to display two arrows pointing in either direction (). The example integral below moves by clicking and dragging the integral center to the right:</p> <p>Before dragging from integral center:  After dragging to the right from integral center: </p>
29	<p>Click Close on the Peak Area/Integration Table when the integration is complete. To zoom back out, right click and click View Entire Spectrum.</p>	<p>The entire spectrum is displayed with integrals:</p> 

Apply a macro

	Action	Result
1	<p>In ProcessIt, click Open Data File icon.</p>  <p>Navigate to "C:\Users\Public\Public Documents\Wiley\KnowItAll\Samples\NMR\Bruker TopSpin\C14H10ClNO3\C14H10ClNO3 H1\fid".</p>	<p>The file is displayed in the main ProcessIt window:</p> 
2	<p>Choose Macro > New/Edit.</p>	<p>The Edit Macro dialog box opens:</p> 

	Action	Result
3	<p>Select H1ProcessSample in the list of Current Macros, then click New.</p>	<p>The Macro Edit Dialog opens:</p>  <p>The steps that are part of the H1ProcessSample macro are displayed on the under Selected functions. The remaining available functions are listed under Available functions.</p>
4	<p>Note: A new macro can be created by using the arrows to remove (<<) or add (>>) items between the two lists of functions in the Macro Edit Dialog.</p> <p>Remove all of the functions from the Selected functions by selecting the item and clicking the left arrow button (<<).</p>	<p>The Selected functions for the Macro Edit Dialog is emptied:</p> 

	Action	Result
5	<p>Move the following functions from Available to Selected using the right arrow button (>>): Default Fourier, Auto Phase, Baseline Correction-Linear, and Auto Integrate.</p> <p>For the Auto Phase function, retain the default parameters (Method = GoodLook). For the Auto Integrate function, change the Noise % value to 1.</p> <p>Note: The macro will execute the functions in the order they are displayed in the Selected section in the Macro Edit Dialog. Items can be reordered using the Up and Down buttons.</p>	<p>The added functions are displayed in order under Selected functions section in the Macro Edit Dialog.</p> 
6	<p>In the Name box in the Macro Edit Dialog, write "New H Macro". Click OK.</p>	<p>The Macro Edit Dialog is closed. "New H Macro" appears on the Edit Macro dialog:</p> 

	Action	Result
7	Click OK on the Edit Macro dialog to close the window. Choose Macro > New H Macro .	<p>The spectrum is processed according to the steps in the applied macro "New H Macro":</p>  <p>The screenshot displays the 'ProcessIt' software window. The title bar reads 'ProcessIt'. The main plot area shows an NMR spectrum for the compound 'C14H10ClNO3 H1'. The x-axis is labeled 'ppm' and ranges from 15 to -1. The y-axis represents intensity, ranging from 0 to 7e+08. The spectrum shows several peaks, with integration values displayed below the baseline: 1.0207, 2.354, 1.646, 1.0000, and 0.329. A prominent peak is visible at approximately 7.2 ppm. The software interface includes a toolbar with various icons for processing and analysis.</p>